

LiftCare

innovative healthcare solutions



The difficult and dangerous process of transferring a patient between bed and stretcher often requires assembling a team of workers who must awkwardly move the patient at great risk to the safety of both patient and caregivers. Current methods frequently result in back injuries, Workcover claims, and patient liability lawsuits—serious problems concerning every healthcare facility.

The Barton® **SAFETTRANS** system is revolutionizing healthcare by providing an easy, safe, and comfortable transfer method with minimum risk to patient and caregiver. This simple, versatile transfer system can help your facility quickly and economically institute a NO-LIFT policy.

The LiftCare Bed Company

Head Office, Showroom & Service Centre, Unit 9, 271 Wells Rd, Chelsea, 3196, Victoria Australia

Tel: +61 3 8773 1111 Fax: +61 3 9773 4008 mail@liftcarebeds.com

www.liftcarebeds.com



Dramatic advances in medical technologies have been made in recent years. One thing yet to be improved, however, is the difficult and awkward process of transferring a patient from bed to stretcher and back again. Patient transfers are a strenuous and dangerous process that frequently results in injury to both patient and caregiver. Often staff must be gathered to accomplish the transfer, and new 'No-Lift' policies are being required.

Tested and proven in hospitals and nursing homes, the *BARTON® SAFETTRANS* reduces both time and caregivers required to transfer a patient, as well as dramatically reducing injuries to both caregivers and patients.

The *BARTON® SAFETTRANS* is design to allow caregivers the ability to transfer a non-ambulatory or limited mobility patient safely and easily. With the *BARTON® SAFETTRANS* a single caregiver can laterally transfer a patient—even one much larger than the caregiver—from almost any hospital bed, procedure/ exam table to stretcher and bed.

Best of all, the *BARTON® SAFETTRANS* allows your facility to implement a NO-LIFT policy quickly and economically. Relatively simple to use, and requiring no stressful or awkward lifting by the caregiver, the *BARTON® SAFETTRANS* uses a simple Set-Clip-and-Crank procedure.

Call us today. The road to safer, more comfortable care—and a NO-LIFT policy for your facility—is just a phone call away.



The System includes:

- 1 Barton® SAFETTRANS bar
- 2 Barton® IV brackets
- 5 Barton® straps with sheet clips and plugs
- 1 Barton® non-slip mattress pad
- 2 Barton® transfer sheets (not shown)



The simplicity of the BARTON® SAFETTRANS allows safe and easy transfer of patients by a single caregiver.



Process showing patient transfer. To transfer back to bed, simply reverses the procedure